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步骤 2: 加入 75 µl of LiDirect™ Extraction Buffer
98°C 孵育 15 min

步骤 3: 加入 125 µl LiDirect™ Stop Solution

步骤 4: set up 25 µl of LiTaq™ Plus PCR reaction mix

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1. 95°C: 15 sec
2. 35 ×
 1. 变性: 5 sec at 98°C
 2. 退火: 10 sec at 55°C
 3. 伸长: 30 sec/kb at 72°C
3. 最终伸长: 5 min at 72°C

产品名称	目录号	规格	价格 (元)
LiDirect™ Lightning Genotyping Kit	M0015-01	80 rxns	¥ 650.00
LiDirect™ Genotyping DNA Extraction Buffer	M0014-04	400 rxns	¥ 1350.00
LiDirect™ Lightning Mouse Genotyping Kit	M0016-01	80 rxns	¥ 650.00
LiDirect™ Lightning Tissue PCR Kit	M0017-01	80 rxns	¥ 650.00
LiDirect™ Lightning Yeast Cell PCR Kit	M0018-01	80 rxns	¥ 650.00

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"The LiDirect™ Lightning Genotyping Kit that I sampled gave clear results and was much more efficient than the methods my lab had been using in the past. With the buffer that is included in the kit, there is no need to make an extraction buffer. Additionally, the quick extraction step (just 15 minutes) really cuts down on time and makes a huge difference. I am very satisfied with the product and plan on using this kit. Thank You!"

----- Icahn SMMS, Emmy Sakakibara, Lab Manager


"LiDirect™ Lightning Genotyping Kit worked very well. But our lab has a lot of previous genotyping kit. After we finish them, I can order your product. Thank you."

----- Brigham and Women's Hospital, Jun Nagai, Postdoc

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
"We were pleasantly surprised at how well the **LiQuant™ Universal Green qPCR Master Mix** worked with our typically problematic samples. I was not familiar with LifeSct prior to this but will definitely keep them in mind for future molecular biology needs and would recommend this product."

----- University of Tennessee, Rebecca Payton, Lab Manager

"I have been using LiTaq™ Plus PCR Mix and **LiQuant™ Universal Green qPCR Master Mix** in my experiments for two years. The PCR mix can impressively amplify the difficult sequences in my hands. Amazingly, it rapidly amplifies the templates at ~ 2kb/min. And the qPCR mix has the higher sensitivity compared with the similar product from other company. I highly recommend these products!!!"


----- UPENN, Janne Liu, Research Associate

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"I have been using LiGreen™ Ultra Nucleic Acid Gel Stain and dNTPs set (100mM each) more than one years, both products are of excellent quality and are cheaper than other companies, I will continue to use them."

----- NIH/NCI, Dr. Ping Yu, Staff Scientist

IntactProtein™ Cell-Tissue Lysis Kit

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"We have tested your IntactProtein™ Cell-Tissue Lysis Kit in liver, adipose tissue, and cultured cells. It worked really well! Extracting proteins from fat is always a headache. Your buffer gave us high protein yields and sharp protein bands. We'd like to switch to your buffer"


----- Columbia University, Dr. Qiang Li, Assistant Professor

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LifeSct LLC, USA

Email: info@lifesct.com
Website: www.lifesct.com

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